



**ADVAMRYL**  
**Poliglecaprone - 25 (Undyed & Violet)**  
 Monofilament Synthetic Absorbable Surgical Suture  
 Supreme handling performance for confidence  
 in subcuticular skin closure



**ADVATRAY™** - the confidence to suture with ease



## Predictable strength for superior support

### Product Description






































Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVA1011	UNDYED	0	70	3/8 CIRCLE		30	-	12
2	ADVA1014	UNDYED	0	90	1/2 CIRCLE		48	-	12
3	ADVA1015	VIOLET	0	70	1/2 CIRCLE		37	-	12
4	ADVA1026	UNDYED	0	70	1/2 CIRCLE		40	HEAVY	12
5	ADVA1027	VIOLET	0	76	1/2 CIRCLE		35	-	12
6	ADVA1029	VIOLET	0	70	1/2 CIRCLE		37	-	12
7	ADVA1030	VIOLET	0	90	1/2 CIRCLE		40	HEAVY	12
8	ADVA1041	VIOLET	0	70	1/2 CIRCLE		40	HEAVY	12
9	ADVA1051	VIOLET	0	75	3/8 CIRCLE		20	-	12
10	ADVA1052	VIOLET	0	70	1/2 CIRCLE		30	-	12
11	ADVA1652	UNDYED	0	70	1/2 CIRCLE		25	-	12
12	ADVA1642	UNDYED	0	70	1/2 CIRCLE		30	-	12
13	ADVA1742	VIOLET	0	70	1/2 CIRCLE		30	-	12
14	ADVA3442	VIOLET	0	70	1/2 CIRCLE		31	LONG TAPER	12
15	ADVA3443	VIOLET	0	70	1/2 CIRCLE		36	LONG TAPER	12
16	ADVA3601	VIOLET	0	70	1/2 CIRCLE		36	-	12
17	ADVA1006	UNDYED	1	90	1/2 CIRCLE		40	-	12
18	ADVA1016	VIOLET	1	70	1/2 CIRCLE		37	-	12
19	ADVA1017	UNDYED	1	90	1/2 CIRCLE		48	-	12
20	ADVA1031	VIOLET	1	70	1/2 CIRCLE		40	HEAVY	12
21	ADVA1042	VIOLET	1	90	1/2 CIRCLE		40	HEAVY	12
22	ADVA1043	VIOLET	1	90	1/2 CIRCLE		50	BLUNT POINT	12
23	ADVA3457	VIOLET	1	70	1/2 CIRCLE		36	LONG TAPER	12
24	ADVA3755	VIOLET	1	70	1/2 CIRCLE		31	LONG TAPER	12
25	ADVA3759U	UNDYED	1	70	1/2 CIRCLE		40	-	12
26	ADVA1005	UNDYED	2-0	70	3/8 CIRCLE		30	-	12
27	ADVA1018	VIOLET	2-0	70	3/8 CIRCLE		24	-	12
28	ADVA1019	VIOLET	2-0	70	1/2 CIRCLE		26	LONG TAPER	12
29	ADVA1020	VIOLET	2-0	70	1/2 CIRCLE		37	-	12
30	ADVA1044	VIOLET	2-0	70	1/2 CIRCLE		35	-	12
31	ADVA1053	VIOLET	2-0	70	1/2 CIRCLE		30	-	12
32	ADVA1075	VIOLET	2-0	75	1/2 CIRCLE		26	LONG TAPER	12

Round Bodied

Oval Round Bodied

Reverse Cutting

Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
33	ADVA1603	UNDYED	2-0	70	1/2 CIRCLE		25	-	12
34	ADVA1602	UNDYED	2-0	70	5/8 CIRCLE		26	-	12
35	ADVA1641	UNDYED	2-0	70	1/2 CIRCLE		30	-	12
36	ADVA1665	UNDYED	2-0	70	1/2 CIRCLE		26	-	12
37	ADVA1666	UNDYED	2-0	70	1/2 CIRCLE		36	-	12
38	ADVA1666S	UNDYED	2-0	45	1/2 CIRCLE		36	-	12
39	ADVA1741	VIOLET	2-0	70	1/2 CIRCLE		30	-	12
40	ADVA3200	VIOLET	2-0	70	1/2 CIRCLE		31	LONG TAPER	12
41	ADVA3327	UNDYD	2-0	70	3/8 CIRCLE		26	-	12
42	ADVA3431	UNDYD	2-0	70	3/8 CIRCLE		26	-	12
43	ADVA3441	VIOLET	2-0	70	1/2 CIRCLE		31	-	12
44	ADVA945P	UNDYED	2-0	90	1/2 CIRCLE		36	-	12
45	ADVA1001	UNDYED	3-0	70	1/2 CIRCLE		20	-	12
46	ADVA1004	UNDYED	3-0	70	3/8 CIRCLE		24	-	12
47	ADVA1008	UNDYED	3-0	90	3/8 CIRCLE		19	-	12
48	ADVA1010	UNDYED	3-0	75	3/8 CIRCLE		26	-	12
49	ADVA1021SA	UNDYED	3-0	70	1/2 CIRCLE		26	-	12
50	ADVA1022	VIOLET	3-0	70	3/8 CIRCLE		24	-	12
51	ADVA1023	VIOLET	3-0	70	1/2 CIRCLE		30	-	12
52	ADVA1028	UNDYED	3-0	70	1/2 CIRCLE		23	-	12
53	ADVA1032	UNDYED	3-0	70	3/8 CIRCLE		24	-	12
54	ADVA1033	VIOLET	3-0	70	1/2 CIRCLE		26	LONG TAPER	12
55	ADVA1035	VIOLET	3-0	70	1/2 CIRCLE		22	LONG TAPER	12
56	ADVA1046	UNDYED	3-0	70	1/2 CIRCLE		26	-	12
57	ADVA1049	VIOLET	3-0	90	3/8 CIRCLE		19	PRIME	12
58	ADVA1050	VIOLET	3-0	75	STRAIGHT		60	PRIME	12
59	ADVA1054	VIOLET	3-0	70	3/8 CIRCLE		20	PRIME	12
60	ADVA1056	UNDYED	3-0	45	3/8 CIRCLE		19	PRIME	12
61	ADVA1064	UNDYED	3-0	45	3/8 CIRCLE		19	-	12
62	ADVA1065	VIOLET	3-0	70	3/8 CIRCLE		19	-	12
63	ADVA1067	UNDYED	3-0	70	3/8 CIRCLE		16	-	12
64	ADVA1070	VIOLET	3-0	75	3/8 CIRCLE		19	-	12
65	ADVA1077	VIOLET	3-0	75	3/8 CIRCLE		26	-	12
66	ADVA1326	UNDYED	3-0	70	3/8 CIRCLE		25	PRIME	12
67	ADVA1326S	UNDYED	3-0	45	3/8 CIRCLE		25	PRIME	12
68	ADVA1326MR	VIOLET	3-0	70	3/8 CIRCLE		25	PRIME	12
69	ADVA1663	UNDYED	3-0	70	1/2 CIRCLE		22	-	12






































 Round Bodied

 Oval Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
70	ADVA1664	UNDYED	3-0	70	1/2 CIRCLE		26	-	12
71	ADVA1737	VIOLET	3-0	70	1/2 CIRCLE		20	-	12
72	ADVA3202	UNDYED	3-0	70	3/8 CIRCLE		19	-	12
73	ADVA3326M	UNDYED	3-0	70	3/8 CIRCLE		26	-	12
74	ADVA3326MS	VIOLET	3-0	45	3/8 CIRCLE		25	-	12
75	ADVA3665	VIOLET	3-0	75	1/2 CIRCLE		31	LONG TAPER	12
76	ADVA1003	UNDYED	4-0	70	3/8 CIRCLE		24	-	12
77	ADVA1007	UNDYED	4-0	75	1/2 CIRCLE		17	-	12
78	ADVA1009	UNDYED	4-0	45	3/8 CIRCLE		19	PRIME	12
79	ADVA1024	VIOLET	4-0	70	3/8 CIRCLE		19	-	12
80	ADVA1025	VIOLET	4-0	70	1/2 CIRCLE		17	-	12
81	ADVA1034	VIOLET	4-0	70	1/2 CIRCLE		22	LONG TAPER	12
82	ADVA1036	UNDYED	4-0	70	3/8 CIRCLE		19	-	12
83	ADVA1038	VIOLET	4-0	75	1/2 CIRCLE		20	LONG TAPER	12
84	ADVA1048	UNDYED	4-0	70	3/8 CIRCLE		19	PRIME	12
85	ADVA1058	VIOLET	4-0	45	3/8 CIRCLE		19	-	12
86	ADVA1066	UNDYED	4-0	70	3/8 CIRCLE		12	-	12
87	ADVA1069	UNDYED	4-0	45	3/8 CIRCLE		13	-	12
88	ADVA1073	VIOLET	4-0	50	3/8 CIRCLE		19	ADVAPASS	12
89	ADVA1076	UNDYED	4-0	45	3/8 CIRCLE		12	-	12
90	ADVA1078	UNDYED	4-0	75	3/8 CIRCLE		20	-	12
91	ADVA1205	UNDYED	4-0	70	3/8 CIRCLE		16	PRIME	12
92	ADVA1205MR	VIOLET	4-0	70	3/8 CIRCLE		16	PRIME	12
93	ADVA1648	UNDYED	4-0	70	3/8 CIRCLE		16	LONG TAPER	12
94	ADVA1748	VIOLET	4-0	70	3/8 CIRCLE		16	LONG TAPER	12
95	ADVA3201	UNDYED	4-0	45	3/8 CIRCLE		19	-	12
96	ADVA3205	UNDYED	4-0	45	3/8 CIRCLE		16	PRIME	12
97	ADVA845G	UNDYED	4-0	45	3/8 CIRCLE		16	PRIME	12
98	ADVA1002	UNDYED	5-0	45	3/8 CIRCLE		19	PRIME	12
99	ADVA1037	VIOLET	5-0	70	1/2 CIRCLE		17	-	12
100	ADVA1039	VIOLET	5-0	75	3/8 CIRCLE		13	LONG TAPER	12
101	ADVA1055	VIOLET	5-0	40	CURVED		20	-	12
102	ADVA1057	UNDYED	5-0	70	3/8 CIRCLE		12	PRIME	12
103	ADVA1059	VIOLET	5-0	45	3/8 CIRCLE		16	-	12
104	ADVA1060	VIOLET	5-0	45	3/8 CIRCLE		16	-	12
105	ADVA1062	UNDYED	5-0	70	3/8 CIRCLE		16	-	12
106	ADVA1068	VIOLET	5-0	70	3/8 CIRCLE		18	-	12

 Round Bodied

 Oval Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
107	ADVA1072	VIOLET	5-0	45	3/8 CIRCLE		12	-	12
108	ADVA1074	VIOLET	5-0	50	3/8 CIRCLE		16	-	12
109	ADVA3203	UNDYED	5-0	45	3/8 CIRCLE		13	PRIME	12
110	ADVA3204	UNDYED	5-0	45	3/8 CIRCLE		16	PRIME	12
111	ADVA3209	UNDYED	5-0	70	3/8 CIRCLE		19	-	12
112	ADVA493G	UNDYED	5-0	45	3/8 CIRCLE		13	PRIME	12
113	ADVA844G	UNDYED	5-0	45	3/8 CIRCLE		16	PRIME	12
114	ADVA1040	VIOLET	6-0	75	1/2 CIRCLE		13	LONG TAPER	12
115	ADVA1061	UNDYED	6-0	70	3/8 CIRCLE		12	PRIME	12
116	ADVA1063	UNDYED	6-0	45	3/8 CIRCLE		13	-	12
117	ADVA1071	UNDYED	6-0	50	3/8 CIRCLE		12	PRIME	12
118	ADVA3214	UNDYED	6-0	45	3/8 CIRCLE		11	PRIME	12
119	ADVA1047	VIOLET	2-0	150	-	-	-	-	12

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
120	ADVA1012	VIOLET	0	75	1/2 CIRCLE			30	12
121	ADVA1021	VIOLET	3-0	70	1/2 CIRCLE & 3/8 CIRCLE			26	12

Round Bodied

Reverse Cutting

Cutting

### Product Description


Composition	Glycolide and ε-Caprolactone
Coating	None
Size	Size 1 to 6-0
Colour	Undyed and Violet
Absorption Profile	Essentially complete up to 90 days
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Trocar, Reverse Cutting, Taper Cut, Blunt Point, Oval Round Bodied & Diamond Point needles

### Tensile Strength

Days of Implantation	Tensile Strength Retention	
	DYED	UNDYED
7 Days (1 to 5-0)	60%	55%
7 Days (6-0)	60%	50%
14 Days (1 to 3-0 & 5-0)	40%	20%
14 Days (4-0 & 6-0)	40%	15%

66% more knot pull tensile strength Vs. recommended limit ensures excellent support throughout the critical wound healing period



### Knot pull tensile strength (Kgf)<sup>1</sup>



REFERENCE \_\_\_\_\_

1 : Data on file with AMS - Bench test done on Poliglecaprone - 25, size 3-0,

### Monofilament Material



### WEEKS



Provides excellent support for a range of critical wound healing period to secure closure in a variety of tissues, procedures and patients.

# ADVA TRAY™ - the confidence to suture with ease

**ADVA TRAY™** a unique dispensing system designed to deliver ease of suture dispensing with minimal memory in the suture strand



## Premium Needles for consistent Superior Performance

**PREMIUM NEEDLES**  
ensure  
**UNSURPASSED SUTURING EXPERIENCE**

- DURABLE COAT**  
RETAINS SHARPNESS...  
PASS AFTER PASS  
Maintains needle sharpness over the duration of procedure, resulting into consistent penetration pass after pass thus less potential for tissue trauma.
- RIBBED BODY**  
FOR SUPERIOR CONTROL  
For superior grasping and stability of needle in needle forceps.
- CHAMFERING**  
MINIMUM SUTURE DETACHMENT  
Diminishes chances for suture fracturing at the swage by rounding off the sharp edges of the barrel in the swage area of the needle.
- ELECTROPOLISHING**  
SMOOTH & BRIGHT  
Makes needle surface free from micro asperities, thus creates ultra sharp edges for lesser tissue drag.

Minimal tissue trauma enables faster recovery



Advanced MedTech Solutions  
P-21-22, 25-26 & 34-35, GIDC, Manjusr, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India  
T : (91) 2667 671302 | [amsltd.com](http://amsltd.com) | [info@amsltd.com](mailto:info@amsltd.com)

Singapore Office : 111, Somerset Road, 07-08, 111 Somerset, Singapore - 238164.  
T: (65) 636910 80 / 89 | E: [info@amsltd.com](mailto:info@amsltd.com) | [amsltd.com](http://amsltd.com)