



Certified Company
Cert No :44 221 127768



Certified Company
Cert No :16 00 B 14030

catalogue

Medical & Hospital Equipment Industries

NURITEK

NURI+EK INDONESIA





CV. Nuri Teknik memulai bisnisnya pada tahun 1991 di atas areal seluas 10.000 m². CV Nuri Teknik memfokuskan usahanya di dalam pembuatan alat-alat Kedokteran dan Rumah Sakit serta Manufaktur Umum.

Pada tahun 1996 CV. Nuri Teknik mendapatkan penghargaan "UPAKARTI" dari presiden RI, selain itu pada tahun 1998 mendapat penghargaan sebagai perusahaan UKM terbaik Jawa Barat dari Gubernur Jawa Barat Bapak R. Nuriana dan pada tahun 2003 mendapatkan penghargaan sebagai Industri Kecil Menengah berprestasi dari Gubernur Jawa Barat Bapak Danny Setiawan. Pada Desember 2006

CV Nuri Teknik terpilih sebagai peraih penghargaan "Dji Sam Soe Award" (tahun 2006 hanya ada 2 pemenang dari 5.112 kandidat). Dengan didukung sarana,teknologi dan tim profesional yang dimiliki, CV. Nuri Teknik terus mengembangkan diri dan telah melakukan ekspor ke Timur Tengah dan Eropa.

CV. Nuri Teknik began business in 1991 on the area 10.000 m². The main products of CV. Nuri Teknik are hospital/medical equipment and also general manufacturer. In 1996 CV. Nuri Teknik got appreciation Upakarti Award from President RI. In 1998 got appreciations the best industry for middle industry from Mr. R. Nuriana The Governor of West Java.

In 2003 got appreciation as the best achievement for middle industry from Mr. Danny Setiawan the Governor of West Java. On December 2006 CV. Nuri Teknik got appreciation "Dji Sam Soe Award" (in year 2006 only two winner from 5.112 candidates). Be supported by high technology instrument and the professional team and now have been exported his product to Middle East and Europe.

PENGHARGAAN YANG PERNAH DIRAIH

"UPAKARTI" dari Presiden Republik Indonesia, pada tahun 1996.

UKM terbaik Jawa Barat dari Gubernur Jawa Barat Bapak R. Nuriana, pada tahun 1998

UKM berprestasi dari Gubernur Jawa Barat Bapak Danny Setiawan, pada Tahun 2003 dan meraih penghargaan "Dji Sam Soe Award" pada Tahun 2006.



UPAKARTI
Appreciation From President
Penghargaan Presiden RI



DJI SAM SOE AWARD

(Tahun 2006 hanya ada 2 pemenang
UKM terbaik dari 5112 kandidat)



Visi Misi

Sejak berdiri CV. Nuri Teknik memiliki visi untuk menjadi Produsen alat kesehatan terbaik di Indonesia. Membuat produk dengan kualitas tinggi dengan harga yang terjangkau dan menyediakan produk di manapun konsumen berada, dalam upaya memenuhi kepuasan pelanggan.

Since the first time CV. Nuri Teknik determined to be the best producer of medical equipments. With high quality and achieved price and give the service where the buyer are. CV. Nuri Teknik always to fulfill the satisfaction of customer.

Produk

Peralatan kedokteran dan rumah sakit yang di hasilkan di CV Nuri Teknik diantaranya : Ranjang Rumah Sakit (dengan berbagai type), Meja Operasi Elektrik, Ginecolog Elektrik, Meja Periksa, Trolley Makan, Lampu Periksa dan lain-lain.

List of product : Hospital bed (with many items of type , Operating table electric, Gynecological electric, Examination table, Food trolley, Examination lamp, and more than 500 items medical and hospital equipment have been product by CV. Nuri Teknik and also we serve the order with design accordance customer request


Dan lebih dari 500 macam alat rumah sakit dan kedokteran yang sudah di produksi CV. Nuri Teknik dan melayani pesanan barang dengan desain sesuai permintaan pemesanan.



Electric Hospital Bed NT 001-DEL4

Specifications

Construction	: Steel pipe
Back Raise	: Adjustable by electric system
Knee Raise	: Adjustable by electric system
H-Low	: Adjustable by electric system
Painting	: Powder Coating
Castor	: 5" two with brake
Dimension	: ± (215L x 98W x 50H) cm
Kubikasi	: ± 2.4m ³



Electric Hospital Bed NT 001-DEL4 (Luxury 2)

Termasuk asesoris : dilengkapi dengan asesoris Side Rail ABS & Roda Central Brake




The above specification are subject to change without prior notice.

Electric ICU Bed NT 001-TRL

Specifications

Construction	: steel pipe
Back Raise	: Adjustable by electric system
Knee Raise	: Adjustable by electric system
Hi-Low	: Adjustable by electric system
Trendelenburg	: Adjustable by electric system
Painting	: Powder coating
Castor	: 5" two with brake
Dimension	: ± (215L x 98W x 60H) cm
Kubikasi	: ± 1.94 m ³



Electric ICU Bed NT 001-TRL (Special Edition)

Termasuk aksesoris : Mattress ABS dilengkapi dengan asesoris Side Rail ABS & asesoris list kaki




The above specification are subject to change without prior notice.

Electric Hospital Bed NT 001-DEL4

Specifications

Construction	: Steel pipe
Back Raise	: Adjustable by electric system
Knee Raise	: Adjustable by electric system
H-Low	: Adjustable by electric system
Painting	: Powder Coating
Castor	: 5" two with brake
Dimension	: ± (215L x 98W x 50H) cm
Kubikasi	: ± 2.4m ³



Electric Hospital Bed NT 001-DEL4 (Grand A1)


Termasuk asesoris : dilengkapi dengan asesoris Side Rail ABS & Roda Central Brake



Electric Hospital Bed NT 001-DEL

Specifications

Construction	: Steel pipe
Back Raise	: Adjustable by electric system
Knee Raise	: Adjustable by electric system
Hi-Low	: Adjustable by electric system
Painting	: Powder Coating
Castor	: 5" two with brake
Dimension	: ± (215L x 98W x 60H) cm
Kubikasi	: ± 2.02 m ³ x 60H) cm



Electric Hospital Bed NT 001-DEL (Prima B1 Plus)

Termasuk asesoris : dilengkapi dengan asesoris Side Rail Aluminum Collapsible




The above specification are subject to change without prior notice.

The above specification are subject to change without prior notice.

Electric Hospital Bed NT 001-DEL4

Specifications

Construction	: Steel pipe
Back Raise	: Adjustable by electric system
Knee Raise	: Adjustable by electric system
H-Low	: Adjustable by electric system
Painting	: Powder Coating
Castor	: 5" two with brake
Dimension	: ± (215L x 98W x 50H) cm
Kubikasi	: ± 2.4m ³



Electric Hospital Bed NT 001-DEL4 (Hilux 1)


Termasuk asesoris : dilengkapi dengan asesoris Side Rail ABS



Electric Hospital Bed NT 001-DEL

Specifications

Construction	: Steel pipe
Back Raise	: Adjustable by electric system
Knee Raise	: Adjustable by electric system
Hi-Low	: Adjustable by electric system
Painting	: Powder Coating
Castor	: 5" two with brake
Dimension	: ± (215L x 98W x 60H) cm
Kubikasi	: ± 2.4 m ³



Electric Hospital Bed NT 001-DEL (Grand A1)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail ABS, aksesoris cover bawah & roda central brake



The above specification are subject to change without prior notice.

The above specification are subject to change without prior notice.

Manual Hospital Bed NT 001-D2

Specifications

Construction : Steel pipe
 Back Raise : Adjustable manually by crank
 Knee Raise : Adjustable manually by crank
 H-Low : Adjustable manually by crank
 Painting : Powder Coating
 Castor : 5" two with brake
 Dimension : ± (215L x 98W x 50H) cm
 Kubikasi : ± 2.02m³



Manual Hospital Bed NT 001-D2

(Prima Manual 1)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminum Collapsible



Manual Hospital Bed NT 001-D2


(Prima Manual 2 Plus)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminium & Telescopic dining table
 Kubikasi: ± 2.02m³

Manual Hospital Bed NT 001-D

Specifications

Construction : Steel pipe
 Back Raise : Adjustable manually by crank
 Knee Raise : Adjustable manually by crank
 Hi-Low : Adjustable manually by crank
 Painting : Powder Coating
 Castor : 5" two with brake
 Dimension : ± (215L x 98W x 60H) cm
 Kubikasi : ± 2.02 m³



Manual Hospital Bed NT 001-D (Prima Manual 1 Plus)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminum Collapsible



Manual Hospital Bed NT 001-C2

Specifications

Construction : Steel pipe
 Back Raise : Adjustable manually by crank
 Knee Raise : Adjustable manually by crank
 Mattress Base : Steel sheet
 Painting : Powder Coating
 Castor : 5" two with brake
 Dimension : ± (215L x 98W x 50H) cm
 Kubikasi : ± 1.21m³

Manual Hospital Bed NT 001-C2

(Prima Manual 1)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminum Collapsible



2
Crank



Manual Hospital Bed NT 001-C2

(Prima Manual 2 Plus)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminum & Telescopic dining table
 Kubikasi: ± 1.21m³

Manual Hospital Bed NT 001-C

Specifications

Construction : Steel pipe
 Back Raise : Adjustable manually by crank
 Knee Raise : Adjustable manually by crank
 Mattress base : steel sheet
 Painting : Powder Coating
 Castor : 5" two with brake
 Dimension : ± (215L x 98W x 60H) cm
 Kubikasi : ± 0.9 m³

Manual Hospital Bed NT 001-C

(Prima Manual 1 Plus)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminum Collapsible



2
Crank



Manual Hospital Bed NT 001-C (Deluxe Manual Plus)


Panel Option : Panel Kayu dilengkapi dengan aksesoris Side Rail stainless steel
 Kubikasi: ± 1.21 m³



Manual Hospital Bed NT 001-B2

Specifications

Construction : Steel pipe
 Back Raise : Adjustable manually by crank
 Mattress Base : Steel sheet
 Painting : Powder Coating
 Castor : 5" two with brake
 Dimension : ± (215L x 98W x 50H) cm
 Kubikasi : ± 0.76m³



Manual Hospital Bed NT 001-B2

(Prima Manual 1)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminum Collapsible



Manual Hospital Bed NT 001-B2

(Prima Manual 2 Plus)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Almunium & Telescopic dining table
 Kubikasi: ± 0.76m³



The above specification are subject to change without prior notice.




Medical & Hospital Equipment Industries



Manual Hospital Bed NT 001-B

Specifications

Construction : Steel pipe
 Back Raise : Adjustable manually by crank
 Mattress base : Steel Sheet
 Painting : Powder Coating
 Castor : 5" two with brake
 Dimension : ± (215L x 98W x 60H) cm
 Kubikasi : ± 0.76 m³



Manual Hospital Bed NT 001-B (Prima Manual 1 Plus)

Termasuk aksesoris : dilengkapi dengan aksesoris Side Rail Aluminum Collapsible



The above specification are subject to change without prior notice.



Medical & Hospital Equipment Industries



Medical & Hospital Equipment Industries

Children Bed NT 014-02CC

Specifications	
Construction	: steel pipe
Head & foot panel	: steel pipe
Back raise	: adjustable manually by crank
Safety side guard	: sliding type, made from steel pipe
Mattress base	: powder coated steel sheet
I.V Stand	: stainless steel pipe
Painting	: powder coating
Castor	: 4" two with brake
Dimension	: ±(175L x 80W x 120H) cm
Kubikasi	: ± 0,68 m ³



Children Bed NT 014-02CCS

Specifications	
Construction	: Stainless steel pipe
Back Raise	: Adjustable manually by crank
Safety side guard	: Sliding type, made from stainless steel pipe
Dimension	: ± (175L x 80W x 120H) cm
Mattress base	: Powder coated steel sheet
I.V Stand	: Stainless steel pipe
Painting	: Powder coating
Castor	: 4", two with brakes
Kubikasi	: ± 0.68m ³



Children Bed NT 014-ICU

(Large) Full Stainless Steel Construction

Specifications	
Construction	: stainless steel pipe
Head & foot panel	: stainless steel pipe & acrylic
Mattress base	: stainless steel sheet
Tredenburg	: adjustable by electric system
Safety side guard	: sliding type, made from stainless steel pipe & acrylic
Dimension	: ± (175L x 80W x 120H) cm
I.V Stand	: stainless steel pipe
Castor	: 4", two with brakes
Kubikasi	: ± 0.68m ³



* Dapat digunakan untuk Children bed dan ICU Children Bed

The above specification are subject to change without prior notice.

Bed Side Cabinet NT 004-1B

Specifications	
Feature	: stainless seel towel hanger, rack shoes, drawer with handle, main cabinet & door with handle
Construction	: steel sheet
Painting	: powder coating
Castor	: 2"
Dimension	: ±(45L x 40W x 85H)cm
Kubikasi	: ±0.32m ³



Bed Side Cabinet NT 004-2A

Specifications	
Top Table	: made from plywood laminated with formica
Feature	: stainless seel towel hanger, rack shoes, drawer with handle, main cabinet & door with handle
Construction	: steel sheet
Painting	: powder coating
Castor	: 2"
Dimension	: ±(45L x 40W x 85H)cm
Kubikasi	: ±0.32m ³



Bed Side Cabinet NT 004-4A

Specifications	
Top Table	: made from plywood laminated with formica
Retractable	: made from plywood laminated with formica
Feature	: stainless seel towel hanger, rack shoes, drawer with handle, main cabinet & door with handle
Construction	: steel sheet
Painting	: powder coating
Castor	: 2"
Dimension	: ±(45L x 40W x 85H)cm
Kubikasi	: ±0.32m ³



Bed Side Cabinet NT 004-4A (Prima)

Specifications	
Top Table	: made from ABS Plastic injection
Retractable	: made from ABS Plastic injection
Feature	: stainless seel towel hanger, rack shoes, drawer with handle, main cabinet & door with handle
Construction	: steel sheet
Painting	: powder coating
Castor	: 2"
Dimension	: ±(45L x 47W x 85H)cm
Kubikasi	: ±0.32m ³



The above specification are subject to change without prior notice.



Over Bed Table NT 005-C

Specifications

Construction : steel pipe
 Adjustable height stand : stainless steel pipe
 Table : plywood laminated with formica
 Painting : powder coating
 Castor : 2"
 Dimension : table : ± (70Lx 40WH) cm
 height : ± (80-120H) cm
 Kubikasi : ± 0,4 m³

Over Bed Table NT 005-B

Specifications

Construction : stainless steel pipe
 Adjustable height stand : stainless steel pipe
 Table : plywood laminated with formica
 Castor : 2"
 Dimension : table : ± (70Lx 40WH) cm
 height : ± (80-120H) cm
 Kubikasi : ± 0,4 m³



Bed Screen NT 022-2A


Specifications

Two part attached together can be fold
 Detachable cloth for cleaning
 Construction : stainless steel pipe
 Castor : 2"
 Dimension : ±(2x90L) x 170H) cm
 Kubikasi : ± 0,38 m³

Hospital Bed Mattress NT 106

Specifications

Made of super foam
 Covered with synthetic leather
 Dimension : ± (200L x 90W x 10H) cm



Children Bed Mattress NT 108-A

Specifications

Made of super foam
 Covered with synthetic leather
 Dimension : ± (10H) cm



Baby Basket Mattress NT 107-A

Specifications

Made of super foam
 Covered with synthetic leather
 Dimension : ± (62L x 32W x 6H) cm



Infusion Stand NT 003-00BSS

Specifications

Feature : adjustable height with locking screw
 Upper stand : stainless steel pipe & rod
 Lower stand : stainless steel pipe, Ø 1 1/4"
 Foot : from plastic star with five castors
 Castor : 2"
 Dimension : height ± (120-220H) cm
 Kubikasi : ± 0,04 m³



Examination Table NT 007-00B

- Specifications**
 Construction : steel pipe
 Back raise : adjustable by manual system
 Mattress : foam with vinyl cover
 Painting : powder coating
 Dimension : ± (200L x 60W x 80H) cm
 Kubikasi : ± 0,47 m³



Examination Table NT 007-00DSS

- Specifications**
 Construction : stainless steel pipe
 Back raise : adjustable by manual system
 Mattress : foam with vinyl cover
 Dimension : ± (200L x 60W x 80H) cm
 Kubikasi : ± 0,47 m³



Foot Step NT 007-F3

- Specifications**
 Construction : stainless steel pipe
 Kubikasi : ± 0,20 m³



Foot Step NT 007-F4

- Specifications**
 Construction : steel pipe
 Painting : powder coating
 Kubikasi : ± 0,20 m³



Examination Stool NT 029-A

- Specifications**
 Construction : stainless steel pipe
 Seat : foam covered with imitation leather
 Foot : star foot type from plastic
 Dimension : ± (Ø35x57-74H) cm
 Kubikasi : ± 0,29 m³



Verlos Bed NT 009-0ASS

- Specifications**
 Head & leg section with detachable lock
 Construction : stainless steel pipe
 Back raise : adjustable by manual system
 Mattress : foam with vinyl cover
 Knee support : adjustable made from stainless steel & aluminum casting
 Handle : stainless steel
 Included : stainless steel pipe I.V stand & stainless steel basin
 Castor : 3" (on leg section)
 Dimension : ± (210L x 80W x 80-120H) cm
 Kubikasi : ± 0,80 m³



Gynecological Examination Table NT 011-00BSS

- Specifications**
 Construction : stainless steel pipe
 Head raise : adjustable by manual system
 Back raise : adjustable by manual system
 Pelvic raise : adjustable by manual system
 Mattress : foam with vinyl cover
 Knee support : adjustable made from stainless steel & aluminum casting
 Included : stainless steel pipe I.V stand & stainless steel basin
 Dimension : ± (195L x 60W x 75H) cm
 Kubikasi : ± 0,63 m³



Gynecological Table /Chair NT 012-00AS

- Specifications**
 Construction : stainless steel pipe
 Back raise : adjustable by manual system
 Leg section : adjustable to horizontal position
 Mattress : foam with vinyl cover
 Included : - stainless steel basin
 - I.V stand stainless steel pipe
 Accessories : knee support
 Dimension : ± (170-180L x 60-70W x 80H) cm
 Kubikasi : ± 0,50 m³



The above specification are subject to change without prior notice.

Examination Lamp NT 023-00B

Steel

Specifications

Feature	: - adjustable position by flexible arm - adjustable height : 130-176H cm - 50 or 100 watt /12V halogen lamp
Upper stand	: steel pipe & flexible arm covered by synthetic rubber
Power consumption	: 50-100 watt /220 VAC
Head lamp & box control	: anodized aluminum extrusion
Foot	: plastic star with five castors
Lower stand	: steel pipe
Painting	: powder coating
Castor	: 2"
Kubikasi	: ± 0,30 m ³



Examination Lamp NT 023-00BSS

Stainless Steel

Specifications

Feature	: - adjustable position by flexible arm - adjustable height : 130-176H cm - 50 or 100 watt /12V halogen lamp
Upper stand	: steel pipe & flexible arm covered by synthetic rubber
Power consumption	: 50-100 watt /220 VAC
Head lamp & box control	: anodized aluminum extrusion
Foot	: plastic star with five castors
Lower stand	: stainless steel pipe
Castor	: 2"
Kubikasi	: ± 0,30 m ³



X-Ray Film Viewer NT 027-00B

Hilux 2, two films

Specifications	
Feature	: screen from acrylic
Lamp	: 2 x 32 watt TL with power switch
Power consumption	: 64 watt /220 VAC
Construction	: steel sheet
Painting	: powder coating
Dimension	: ± (75L x 12W x 55H) cm
Kubikasi	: ± 0,16 m ³

Resusitasi Table NT 007-EB

Specifications

Completed with	: 2 drawers and lockers
Construction	: steel sheet
Painting	: powder coating
Dimension	: ± (100L x 60W x 80H) cm
Kubikasi	: ± 0,77 m ³



Emergency Stretcher NT 015-01 BSS

Specifications

Construction	: stainless steel pipe
Safety side guard	: stainless steel pipe
Back raise	: adjustable by manual system
Completed	: with rack & oxygen rack
I.V stand	: stainless steel pipe & rod
Mattress	: foam with vinyl cover
Castor	: 5" two with brakes
Dimension	: ± (200L x 60W x 85H) cm
Kubikasi	: ± 0,96 m ³



Stretcher NT 015-00BSS

Specifications

Construction	: stainless steel pipe
Back raise	: adjustable by manual system
I.V stand	: stainless steel pipe & rod
Mattress	: foam with vinyl cover
Castor	: 5" two with brakes
Dimension	: ± (200L x 60W x 85H) cm
Kubikasi	: ± 0,57 m ³



Stretcher NT 015-00B

Specifications

Construction	: steel sheet
Back Raise	: adjustable by manual system
I.V stand	: stainless steel pipe & rod
Mattress	: foam with vinyl cover
Painting	: powder coating
Castor	: 5" two with brakes
Dimensions	: ± (200L x 60W x 85H) cm
Kubikasi	: ± 0,57M ³

Instrument Trolley NT 016-B

Stainless Steel

Specifications
 Two shelves
 Construction : stainless steel pipe
 Shelf : stainless steel sheet
 Guard list : stainless steel pipe
 Handle : stainless steel pipe
 Painting : powder coating
 Castor : 2"
 Dimension : ± (60L x 40W x 80H) cm
 Kubikasi : ± 0,42 m³



Instrument Trolley NT 016-J

Stainless Steel

Specifications
 Two shelves
 Knock down system
 Construction : stainless steel pipe
 Shelf : stainless steel sheet
 Castor : 2"
 Dimension : ± (60L x 40W x 80H) cm
 Kubikasi : ± 0,12 m³



Dressing Trolley NT 018-02ASS

Stainless Steel

Specifications
 Stainless steel basin
 Plastic bucket
 Construction : stainless steel pipe
 Upper shelf : stainless steel sheet, consist of two drawers with stainless steel pipe as list
 Lower shelf : stainless steel pipe as list
 Drawer : stainless steel sheet
 Handle : stainless steel pipe
 Castor : 2"
 Dimension : ± (75L x 50W x 85H) cm
 Kubikasi : ± 0,51 m³



The color of plastic bucket depend to our stock

Mayo Table NT 006-B

Specifications
 Construction : stainless steel pipe
 Adjustable height stand : stainless steel pipe
 Tray : stainless steel sheet
 Castor : 2"
 Dimension : table : ± (60L x 40H) cm
 height : ± (80- 120H) cm
 Kubikasi : ± 0,4 m³



Emergency Trolley NT 019-01 BSS

Specifications
 Three drawers
 Two locker
 Adjustable side table
 Oxygen rack stainless steel rod
 I.V stand stainless steel pipe & rod
 Construction : stainless steel sheet
 Top table : stainless steel sheet
 Castor : 4" two with brakes
 Dimension : ± (65L x 45W x110H) cm
 Kubikasi : ± 0,66 m³



Bowl Stand NT 017-2A

Specifications
 Star foot type from plastic
 Construction : stainless steel pipe
 Basin : two pieces, stainless steel
 Towel hanger : stainless steel rod
 Castor : 2"
 Dimension : ± (75L x 76H) cm
 Kubikasi : ± 0,25 m³



Medication Cart/Trolley

NT 031-A

Specifications

One drawer	
Push handle	
Two doors with lock and handle	
Construction	: stainless steel sheet
Medicine drawer	: 32 pieces
Castor	: 4" two with brakes
Dimensions	: ± (70L x 40W x 100H) cm
Kubikasi	: ± 0,59 m ³



Ultra Violet Sterilization

NT 101-00A

Deluxe, stainless steel

Specifications

Features	: - the most simply and effective for sterilize operation room - UV light sterilizer room from all bacterial
Completed	: with timer, hour meter & castor
Construction	: stainless steel pipe & sheet
Power supply	: 220 V AC/50HZ
Power consumption	: 180 Watt
UV lamp	: 30 Watt x 4 pieces, 10 Watt x 1 pieces
Lamp durability	: 3.000 hours
Dimensions	: ± (45L x 45W x 130H) cm
Kubikasi	: ± 0,51 m ³



Phototherapy/Blue Light

NT 032-00A

Specification

Construction	: square steel pipe & steel sheet
Power supply	: 220 V AC,50Hz
Power consumption	: 150 Watt
Lamp	: TL blue 20 Watt x 6 Pcs
Hour meter	: 0~9999 hours
Painting	: powder coating
Castor	: 4" two with brakes
Dimensions	: - head : ± (70L x 70W x 10H) cm - base : ± (70L x 70W x 100~155H) cm
Kubikasi	: ± 0,66 m ³



Mortuary Cabinet

NT 112-2J

Specifications

Unit	: 1 (one) unit capacity for 2 corpses
Temperature	: (+) 2°C ~ (+) 8°C
Isolation	: polyurethane foam thickness 7,5 cm
Body	: stainless steel sheet aisi 304
Doors	: 2 (two) unit swing door
Machine	: 1 (one) unit
Evaporator	: 1 (unit) blower system
Refrigerant	: R-22
Equipments	: thermostat, expansion valve, panel electric, filter dryer, thermometer, 2 unit corpse stretcher, etc.
Dimension	: ± (280L x 90W x 190H) cm



Automatic Scrub Station NT 111-B

Specifications

Construction	: stainless steel sheet
Power consumption	: 300 watt
Power supply	: 220V/50HZ
Minimum water pressure	: 2 kg/cm ²
Adjustable flow volume	: max. 15 ltr/minutes
Sterilizing method	: sterilization by UV lamp
Filtering	: 05, 02 for pre filtering and 01 micron for main filter
Dimension	: ± (140L x 60W x 190H) cm
Kubikasi	: ± 1,55 m ³



Automatic Scrub Station NT 111-C

Specifications

Construction	: stainless steel sheet
Power consumption	: 300 watt
Power supply	: 220V/50HZ
Minimum water pressure	: 2 kg/cm ²
Adjustable flow volume	: max. 15 ltr/minutes
Sterilizing method	: sterilization by UV lamp
Filtering	: 05, 02 for pre filtering and 01 micron for main filter
Dimension	: ± (200L x 60W x 190H) cm
Kubikasi	: ± 2,30 m ³

Instrument Cabinet NT 025-1B

Specifications

Single door with lock and handle	
Three shelves	
Construction	: steel sheet
Painting	: powder coating,
Dimensions	: ± (65L x 45 W x 170H) cm
Kubikasi	: ± 0,79 m ³



Instrument Cabinet NT 025-2B

Specifications

Double door with lock and handle	
Three shelves	
Construction	: steel sheet
Painting	: powder coating,
Dimensions	: ± (120L x 45 W x 170H) cm
Kubikasi	: ± 1,50 m ³

Medicine Cabinet NT 025-1D (Deluxe)

Specifications

Single door with lock and handle	
Three shelves	
One lockers with lock and handle	
Construction	: steel sheet
Painting	: powder coating
Dimensions	: ± (65L x 45W x 170H) cm
Kubikasi	: ± 0.79M ³



The above specification are subject to change without prior notice.

The above specification are subject to change without prior notice.

Baby Basket

NT 013-00BSS

Stainless Steel & Acrylic

Specifications

Basket can be moved and tilt in trendelenburg position

- Construction : stainless steel pipe
- Basket : stainless steel rod
with acrylic on every side
- Castor : 2"
- Dimension : Trolley : $\pm(70L \times 41W \times 72H)$ cm
Basket : $\pm(64L \times 35W \times 35H)$ cm
- Kubikasi : $\pm 0,32 \text{ m}^3$



Baby Basket

NT 013-00A

Steel

Specifications

Basket can be moved and tilt in trendelenburg position

- Construction : steel pipe
- Basket : steel rod
with steel wiremesh on every sides
- Castor : 2"
- Painting : powder coating
- Dimension : Trolley : $\pm(70L \times 41W \times 72H)$ cm
Basket : $\pm(64L \times 35W \times 35H)$ cm
- Kubikasi : $\pm 0,32 \text{ m}^3$



Laundry Trolley

NT 030-A

Specifications

- Cloth
- Push handle : steel pipe
- Construction : angle steel
- Painting : powder coating
- Castor : 4" two with brakes
- Dimensions : $\pm(90L \times 60W \times 75H)$ cm
- Kubikasi : $\pm 0,78 \text{ m}^3$



Narcotic Cabinet

NT 025-NCA

Specifications

- Construction : steel sheet
- Feature : - double doors, made from steel plate thickness 3mm
each door with lock (inner and outdoor lock)
- ready installed in the wall
- Painting : Powder coating
- Dimension : $\pm(60L \times 30W \times 120H)$ cm
- Kubikasi : $\pm 0.36 \text{ M}^3$



Narcotic Cabinet

NT 025-NCB

Specifications

- Construction : steel sheet, thickness 1.5 mm
- Feature : - double doors
- ready installed in the wall
- Painting : Powder coating
- Dimension : $\pm(70L \times 40W \times 80H)$ cm
- Kubikasi : $\pm 0.36 \text{ M}^3$

Narcotic Cabinet

NT 025-NCC

Specifications

- Construction : steel sheet, thickness 2 mm
- Feature : - double doors
- ready installed in the wall
- Painting : Powder coating
- Dimension : $\pm(90L \times 60W \times 120H)$ cm
- Kubikasi : $\pm 0.91 \text{ M}^3$



Gynecological Chair NT 036-00A2SS

Specifications

Feature :

- Easy adjustable in height by remote control
- Easy adjustable in (reverse) trendelenburg position and back raise position by remote control
- Hi-low (reverse) trendelenburg and back raise adjustable by electric system

Chassis construction : Steel

Body cover : Stainless Steel sheet

Mattress : Foam with vinyl cover

Painting : Powder coating

Accessories :

- knee support : Adjustable, made from stainless steel and aluminum casting
- hand holder : Adjustable, made from stainless steel
- stainless steel basin (square)
- Power supply : 220V AC/50Hz
- Dimension : ± (190L x 60W x min 70 - max. 95H) cm



Operating Table NT 033-00CELSS

Specifications

Feature adjustable head support, backrest and leg support

All the body is divided to be four sections : head support, back, center and leg part

Back raise,tilting, hi-low, trendelenburg / reverse trendelenburg by electric system

Construction : stainless steel sheet

Mattress : foam with vinyl cover

Dimension : ±(200L x 60W x ±min70H ~ ±max 100H) cm



Hospital Bed Multipurpose Bed/Bed Partus NT 037-B

Specifications

Adjustable leg section

Can be transform in 2 mode, hospital bed mode & gynecological chair mode

Construction : steel pipe

Mattress base : acrylic or steel sheet

Back raise : adjustable by electric system

Hi-low : adjustable by electric system

Trendelenburg : adjustable manual, use double gas spring system

Head & foot panel : ABS & detachable

Safety side guard: ABS, with smooth mechanism used gas spring

Basin : stainless steel

Knee support : adjustable, made form stainless steel & alumunium casting

Foot step : can be attached or detached

Include : - hand grip

- stirrup

- rack for accessories hanger

Mattress : foam with vinyl cover

Painting : powder coating

Castor : 5" two with brakes

Dimension : ± 165-220L x 92W x 55-85H cm



Gynecological Chair Mode



Rack Accessories



Patient Transport/Transfer Bed Hydraulic NT 015-TBH

Specifications

Construction : stainless steel pipe, steel pipe,
stainless steel sheet & steel sheet

- Head & Foot panel
- Side Rail
- Oxygen Rack
- 4 pieces Castors
- Use bumper on each sides

Back raise : adjustable by gas spring
Hi-Low : adjustable by hydraulic
Trendelenburg : adjustable by manually
Painting : powder coating
Dimension : ±(205L x 68W x 53 - 86H)cm




Produk Baru Tayang di E-Katalog

Patient Transport NT 015-C

Specifications

Construction : stainless steel pipe, steel pipe & steel sheet
Head & foot panel : stainless steel pipe
Back raise : adjustable by electric system
Hi-low : adjustable by electric system
Mattress : foam with vinyl cover
Painting : powder coating
Oxygen rack : stainless steel rod
I.V stand : stainless steel pipe
Castor : 5" two with brakes
Dimension : ±(200L x 60W x 55-85H) cm



Gynecological Examination Table NT 011-00FEL

Specifications

Stainless steel basin
Construction : steel pipe
Pelvic rest : adjustable by electric system
Head rest : adjustable by electric system
Back rest : adjustable by electric system
Hi-low : adjustable by electric system
Knee support : adjustable, made from stainless steel
& aluminium casting
I.V. stand : stainless steel pipe & rod with detachable
Mattress : foam with vinyl cover
Painting : powder coating
Dimension : ± 195L x 60W x 50 - 80H cm



Examination Table NT 007-00GS

Specifications

Construction : steel pipe
Back raise : adjustable by electric system
Mattress : foam with vinyl cover
Hi-low : adjustable by electric system
Painting : powder coating
Dimension : ± (200W x 60W x 55-85H) cm



Examination Table NT 007-F

Specifications

Adjustable back raise
Construction : steel pipe
Mattress : foam with vinyl cover
Painting : powder coating
Dimension : ± (200L x 60W x 80H) cm

